

.....Experiment Summary.....
A series of experiments was conducted to measure the potential of a quinone solution at various pH levels, specifically from pH 7 to pH 9. The procedure involved adding 6 mL of buffer solution of the desired pH, 3 mL of NaCl solution, 15 mL of water, and 6 mL of quinone solution to a beaker. A CV scan was then run to measure the potential. The contents of the beaker were then transferred to a waste beaker. The observations showed that as the pH increased, the potential became more negative. The potential at pH 3.87 was -0.388 eV, at pH 4.79 was -0.443 eV, at pH 6.02 was -0.505 eV, at pH 6.95 was -0.576 eV, at pH 7.85 was -0.631 eV, and at pH 9.13 was -0.674 eV.

After performing the experiments, these are the results:

The estimate for pKa1 is 8.096.
The estimate for pKa2 is 12.380.
The estimate for slope is -60.958.

