

RObust biOmarker Feature Selection (ROOFS)

Step 1: Pipeline configuration

Required

> Dataset

> Feature selection methods

> Configuration

> Outputs

Optional

> Fixed signature

Step 2: Run the pipeline

Run

▼ Feature selection methods

The default selected methods are the methods recommended in the file README_methods.md. Some methods require to set a fixed/min/max number of features to be selected. Others need specification of a number of clusters. If selecting such methods, the value of these hyperparameters need to be set in the Configuration section.

no specification

fixed set size

minimal set size

maximum set size

minimum and maximum set sizes

number of clusters

> Core methods

▼ Extended methods

FS methods available in the full version of the package, which requires a more advanced installation process.

Embedded

☐ Adaptive Lasso

☐ Exclusive Lasso

☐ Stochastic Gates

☐ Hilbert Schmidt Independence Criterion Lasso (HSIC Lasso)

☐ LassoNet

☐ Relaxed Lasso

☐ Repeated Elastic Net Technique

☐ Sparse group Lasso, correlation distance metric

☐ Sparse group Lasso, Minkowski distance metric

☐ Stability Selection